

Memo

To: AML Reference Laboratory Clients
From: Allina Medical Laboratories, Reference Laboratory Services
Date: March 17, 2006
Re: Amylase, Timed Urine Panels Discontinued

Effective immediately, the following panels have been made non-orderable:

<i>33/AM2</i>	<i>Amylase, 2 Hour Urine</i>
<i>36/AMU</i>	<i>Amylase, 24 Hour Urine</i>

For any timed urine amylase request, the following panel should be ordered:

Test Name:	AMYLASE, TIMED URINE
Test Number:	29/ATU
Collect:	10 ml Timed Urine
Container:	Plastic Urine Container
Processing:	Submit Aliquot or Entire Specimen
Transport/Stability:	Refrigerated
Performing Lab:	AML
Days Set Up:	M - Su
Expected TAT:	1 Day
Ref. Ranges:	See Report
Collection/ Processing Details:	Record Volume and number of hours of collection on request and aliquot.
Method:	Enzymatic
CPT Codes:	82150