

Memo

To: AML Reference Laboratory Clients
From: Allina Medical Laboratories, Reference Laboratory Services
Date: May 23, 2007
Re: Estrone Specimen Requirement Change Notification

AML has received notification from Quest Diagnostics that **effective immediately**, the required specimen for Estrone (4763/ERO) is 1ml of serum from a plain red top tube.

The AML on line Collection Manual, located at www.allina.com/medicallaboratories has been updated to reflect this change

Test Name:	ESTRONE
Test Number:	4763/ERO
Collect:	1 ml Serum - Plain Red
Alternate Collect:	NONE
Container:	Plastic Screw Cap Transport Vial
Processing:	Spin and Separate
Transport/Stability:	Refrigerated
Alternate Names:	E1
Performing Lab:	Quest Diagnostics
Days Set Up:	T - Sa
Expected TAT:	3 - 6 Days
Ref. Ranges:	See Report
Collection/	
Processing Details:	Rejected if on SST, Gross Hemolysis or Lipemic
Method:	LC/MS/MS
CPT Codes:	82679-90
Date Created:	11/30/2001
Revised Date:	05/23/2007