



***PINWORM PREP  
6603/PIN***

***Patient Collection Instructions:***

**Note:** The optimal time to collect specimens for pinworm examination is upon arising in the morning, before bathing or defecating.

- Remove the cap with the clear plastic paddle from the collection container. Do not touch the paddle.
- Press the sticky surface of the paddle (marked "sticky side") against the perianal skin folds with moderate pressure.
- Return the paddle/cap to the container.
- Print the patient's name and date of birth on the container label.
- Return the container to your doctor's office/clinic at room temperature.